

# MDW-160M Mechanical Person Scale



## Technical Specifications

Maximum capacity	350lb/160kg
Readability	4oz/100g
Repeatability (S.D.)	0.4lb / 200g
Linearity ( $\pm$ )	0.4lb / 200g
Platform size	14.8x10.8" / 375x275mm

## General Specifications

Model	MDW-160M
Units of measurement	Kilogram (kg), Pound (lb)
Housing	Metal
Overall dimensions (wxdxh)	20.9x10.8x58.5" / 530x275x1485mm



Adam Equipment - offices worldwide

[www.adamequipment.com](http://www.adamequipment.com)