




MDW-160M Mechanical Person Scale



Speed, Performance, Value - Perfect Balance

- Height rod as standard (measures up to 60-212 cm)
- Robust metal housing
- Anti-slip platform
- Adjustable zero point
- Wheels to allow easy movement
- 2 weighing units (kg, lb)
- Weighing
- Measuring height




Technical Specifications

Maximum capacity	350lb/160kg
Readability	4oz/100g
Repeatability (S.D.)	0.4lb / 200g
Linearity (\pm)	0.4lb / 200g
Platform size	14.8x10.8" / 375x275mm

General Specifications

Model	MDW-160M
Units of measurement	Kilogram (kg), Pound (lb)
Housing	Metal
Overall dimensions (wxdxh)	20.9x10.8x58.5" / 530x275x1485mm



Adam Equipment - offices worldwide

www.adamequipment.com