





Speed, Performance, Value - Perfect Balance

- Height rod as standard (measures up to 71-210 cm)
- Robust metal housing
- Anti-slip platform
- Large, clear backlit LCD display
- Colour coded, sealed keypad
- External calibration
- Hold function
- Automatic backlight
- Full range tare
- AC Adapter
- Overload protection
- Rechargeable battery / mains adapter supplied as standard
- 2 weighing units (kg, lb)
- Precision load cell technology
- Simple user-friendly operation
- Weighing
- Measuring height



Technical Specifications

Maximum capacity	550lb/250kg
Readability	0.2lb/100g
Repeatability (S.D.)	0.4lb / 0.2kg
Linearity (±)	0.4lb / 200g
Platform size	14.6x10.8" / 372x275mm

General Specifications

Model	MDW-250L
Tare functions	Yes
Units of measurement	Kilogram (kg), Pound (lb)
Applications	Weighing and measuring height
Functions	Hold and Tare functions
Calibration	External
Display	LCD
Power supply	Internal re-chargeable 6V battery / 220 VAC through adapter
Operating temperature	41°F - 104°F
Height range	71cm - 210cm
Division of measurement	0.5cm
Other features	Automatic backlight, overload protection
Housing	Metal, Indicator: ABS Plastic
Overall dimensions (wxdxh)	21.1x10.8x37" / 535x275x940mm



Adam Equipment - offices worldwide

www.adamequipment.com