

Technical Specifications


Model	TBB 2610T	TBB 2610S	TBB 610S
Maximum capacity	2610g	2610g	610g
Readability	0.1g	0.1g	0.1g
Tare range	+225g	-	-
Repeatability (S.D.)	0.2g	0.2g	0.2g
Linearity (\pm)	0.2g	0.2g	0.2g
Pan size	6" \varnothing / 152mm \varnothing		
Net weight	5.5lb		3.8lb

General Specifications

Units of measurement	Gram (g)
Calibration	External
Display	Graduated Beam
Housing	Aluminium
Overall dimensions (wxdxh)	19.9x4.3x6.3" / 505x110x160mm

Key Accessories

- Security Lock and Cable
- Hard Carry Case with Lock



Adam Equipment - offices worldwide

www.adamequipment.com