

Technical Specifications


| Model | TBB 2610T | TBB 2610S | TBB 610S |
|----------------------|----------------|-----------|----------|
| Maximum capacity | 2610g | 2610g | 610g |
| Readability | 0.1g | 0.1g | 0.1g |
| Tare range | +225g | - | - |
| Repeatability (S.D.) | 0.2g | 0.2g | 0.2g |
| Linearity (±) | 0.2g | 0.2g | 0.2g |
| Pan size | 6" Ø / 152mm Ø | | |
| Net weight | 5.5lb | | 3.8lb |

General Specifications

| | |
|----------------------------|--|
| Units of measurement | Gram (g) |
| Applications | Weighing, Net weight / tare, Below balance weighing, Density of liquids or solids, Filling |
| Calibration | External |
| Display | Graduated Beam |
| Housing | Aluminium |
| Overall dimensions (wxdxh) | 19.9x4.3x6.3" / 505x110x160mm |

Key Accessories

- Security Lock and Cable
- Hard Carry Case with Lock



Adam Equipment - offices worldwide

www.adamequipment.com