

Technical Specifications

General Specifications

Model	AE401
Tare functions	Full
Interface	RS-232 Uni-directional
Operating temperature	32°F - 104°F
Power supply	12VDC 800ma adaptor
Display	7 digit LCD with backlight
Housing	Stainless Steel
Overall dimensions (wxdxh)	9.4x5.7x3.7" / 240x145x95mm

Key Accessories

- Adam Printer
- Thermal Paper for Adam Printer (5 Pack)



Adam Equipment - offices worldwide

www.adamequipment.com