

MDW-250L Person Scale

Technical Specifications

Maximum capacity	550lb/250kg
Readability	0.2lb/100g
Repeatability (S.D.)	0.4lb / 0.2kg
Linearity (\pm)	0.4lb / 200g
Platform size	14.6x10.8" / 372x275mm

General Specifications

Model	MDW-250L
Tare functions	Yes
Units of measurement	Kilogram (kg), Pound (lb)
Applications	Weighing and measuring height
Functions	Hold and Tare functions
Calibration	External
Display	LCD
Power supply	Internal re-chargeable 6V battery / 220 VAC through adapter
Operating temperature	41°F - 104°F
Height range	71cm - 210cm
Division of measurement	0.5cm
Other features	Automatic backlight, overload protection
Housing	Metal, Indicator: ABS Plastic
Overall dimensions (wxdxh)	21.1x10.8x37" / 535x275x940mm



Adam Equipment - offices worldwide

www.adamequipment.com