



## Speed, Performance, Value - The Right Balance

- Magnetic damping for faster stability
- Rigid construction for durability
- Zero adjust facility to enable zero reading
- Notched and tiered graduated beams with centre reading
- Large stainless steel weighing pan
- Built-in platform for density determination
- Weighing
- Below balance weighing




## Technical Specifications

Maximum capacity	311g
Readability	0.01g
Repeatability (S.D.)	0.02g
Linearity ( $\pm$ )	0.02g
Pan size	100mm $\varnothing$

## General Specifications

Model	AQ 311S
Units of measurement	Gram (g)
Functions	Weighing
Calibration	External
Display	Graduated Beam
Other features	Magnetic damping, Density platform with adjustable height
Housing	Aluminium
Overall dimensions (wxdxh)	370x300x137mm
Net weight	1.3kg



Adam Equipment - offices worldwide

[www.adamequipment.com](http://www.adamequipment.com)