

Technical Specifications


Model	BFW 75	BFW 150	BFW 300
Maximum capacity	75kg	150kg	300kg
Readability	5g	10g	20g
Repeatability (S.D.)	5g	10g	20g
Linearity (\pm)	10g	20g	40g
Platform size	560x470mm		

General Specifications

Tare functions	Full
Units of measurement	Gram (g), Kilogram (kg), Pound (lb), Ounce (Oz)
Applications	Weighing, check weighing
Typical stabilisation time	3 - 5 Sec
Interface	RS-232 bi-directional
Calibration	Automatic calibration
Display	Large LCD with backlight
Power supply	Internal re-chargeable battery or mains power
Operating temperature	0°C to 40°C
Housing	Indicator: ABS Plastic, Platform: Mild Steel base with stainless steel pan
Overall dimensions (wxdxh)	470x685x950mm
Net weight	20kg

Key Accessories

- Adam Printer
- AdamDU - Data Collection Program
- Thermal Paper Roll for Adam Printer
- In-use Cover
- RS-232 Cable



Adam Equipment - offices worldwide

www.adamequipment.com