



The Adam Printer is a compact thermal printer which is ideally suited for use with the Adam range of balances and scales.

The printer is controlled directly from the balance and is supplied with an RS-232 connector which enables it to be linked directly to the balances. As a general purpose printer it can be used with any RS-232 interface.

The Adam Printer will also record full calibration details in compliance with GLP requirements when used with the Adam balances. When used with the PMB moisture analyser, the printer will record the % loss in weight at each preselected interval, and at the conclusion of the test prints a record of the test parameters.

Power Supply	8.5 to 14VDC or 7 to 10VAC from module 15watts when printing 3watts when idle
Mains Supply	As required 230VAC 50/60hZ
Environment	35°C to 52°C 10% to 80% humidity, non condensing
Overall Dimensions	160mm x 130mm x 45mm
Net Weight	0.5kg

## Accessories

- Item No. 9192 - Thermal Replacement Paper (5-Pack)