


## Technical Specifications

Maximum capacity	311g
Readability	0,01g
Repeatability (S.D.)	0,02g
Linearity ( $\pm$ )	0,02g
Pan size	100mm $\varnothing$

## General Specifications

Model	AQ 311S
Units of measurement	Gram (g)
Functions	Weighing
Calibration	External
Display	Graduated Beam
Other features	Magnetic damping, Density platform with adjustable height
Housing	Alumunium
Overall dimensions (wxdxh)	370x300x137mm
Net weight	1,3kg



Adam Equipment - offices worldwide

[www.adamequipment.com](http://www.adamequipment.com)