

Technical Specifications


| Model | CBW 3 | CBW 6H | CBW 6 | CBW 15 | CBW 30 | CBW 45 |
|----------------------|-----------|---------|---------|---------|---------|---------|
| Maximum capacity | 3000g | 6000g | 6000g | 15kg | 30kg | 45kg |
| Readability | 0.1g | 0.1g | 0.2g | 0.5g | 1g | 2g |
| Resolution | 1:30000 | 1:60000 | 1:30000 | 1:30000 | 1:30000 | 1:22500 |
| Tare range | -3kg | -6kg | -6kg | -10kg | -30kg | -45kg |
| Repeatability (S.D.) | 0.1g | 0.1g | 0.2g | 0.5g | 1g | 2g |
| Linearity (±) | 0.2g | 0.2g | 0.4g | 1g | 2g | 4g |
| Pan size | 225x275mm | | | | | |

General Specifications

| | |
|----------------------------|--|
| Units of measurement | Gram (g), Kilogram (kg) |
| Functions | Weighing, basic parts counting, % weighing and checkweighing |
| Typical stabilisation time | 2 Sec |
| Calibration | Automatic External |
| Display | 1 x 6 digits LCD digital display |
| Power supply | 9 VDC, 800 mA from external adapter |
| Operating temperature | 0°C - 40°C |
| Housing | ABS plastic, stainless steel platform |
| Overall dimensions (wxdxh) | 315x355x110mm |
| Net weight | 4.1kg |

Key Accessories

- In-use Cover
- Adam Printer
- In-Use Cover pack of 10
- In-Use Cover pack of 5
- Weigh Below Hanger (factory fitted)
- Thermal Paper Roll for Adam Printer



Adam Equipment - offices worldwide

www.adamequipment.com