



## Speed, Performance, Value - The Right Balance

- Large LCD display
- Large stainless steel pan
- Sealed keypad
- Tare function
- 5 weighing units (kg, g, lb, oz, lbs:oz)
- Overload protection
- Non-slip adjustable levelling feet
- Auto sleep / power down function to save battery life
- Colour coded, sealed keypad
- Easy access to the rechargeable battery
- External calibration
- Full range tare
- Large easy to read LCD display
- Low battery indication
- AC Adapter
- Precision load cell technology
- Rechargeable battery up to 150 hours operation
- Simple user-friendly operation
- Splashproof to protect from accidental spills
- Zero Tracking
- Weighing



## Technical Specifications

Model	QBW 1500	QBW 3000	QBW 6000	QBW 15	QBW 30
Maximum capacity	1500g	3000g	6000g	15kg	30kg
Readability	0.5g	1g	2g	5g	10g
Resolution	1:3000	1:3000	1:3000	1:3000	1:3000
Repeatability (S.D.)	0.5g	1g	2g	5g	10g
Linearity ( $\pm$ )	0.5g	1g	2g	5g	10g
Pan size	225x275mm				

## General Specifications

Units of measurement	Gram (g), Kilogram (kg), Pound (lb), Ounce (Oz), Pound-Ounce (lb-oz)
Applications	General purpose weighing scale
Typical stabilisation time	2 Sec
Calibration	Manual external
Display	Large LCD display with 25mm/1" digits
Power supply	Rechargeable battery / mains adapter (adapter supplied as standard)
Operating temperature	0°C - 40°C
Housing	ABS Plastic, stainless steel platform
Overall dimensions (wxdxh)	315x355x110mm
Net weight	4.1kg

## Key Accessories

- Weigh Below Hanger (factory fitted)
- In-use Cover
- In-Use Cover pack of 10
- Hard Carry Case with Lock
- In-Use Cover pack of 5



Adam Equipment - offices worldwide

[www.adamequipment.com](http://www.adamequipment.com)