

# CBC M Bench Counting Scales (EC Approved)



## Speed, Performance, Value - The Right Balance

- EC type approval
- Non-slip adjustable levelling feet
- Large backlit LCD display
- Overload protection
- Easy to use keypad
- Rechargeable battery up to 90 hours operation
- Preset tare function
- Large stainless steel pan
- Memory accumulation
- Auto sleep / power down function to save battery life
- Colour coded, sealed keypad
- Date and time
- Easy access to the rechargeable battery
- Full colour coded keypad with numeric entry
- Full range tare
- RS-232 bi-directional interface
- Simple user-friendly operation
- Splashproof to protect from accidental spills
- Adjustable filters
- Zero Tracking
- Internal counting resolution 1:400,000
- Low battery indication
- AC Adapter
- Precision load cell technology
- Pre-set counting with alarm
- Weight accumulation
- Parts counting
- Check counting
- Count accumulation



## Technical Specifications

Model	CBC 3M	CBC 6M	CBC 15M	CBC 30M
Maximum capacity	3000g	6000g	15kg	30kg
Readability	1g	2g	5g	10g
Resolution	1:3000	1:3000	1:3000	1:3000
Tare range	-3kg	-6kg	-10kg	-30kg
Repeatability (S.D.)	1g	2g	5g	10g
Linearity (±)	1g	2g	5g	10g
Pan size	225x275mm			

## General Specifications

Internal counting resolution	1:300,000
Units of measurement	Kilogram (kg)
Typical stabilisation time	2 Sec
Interface	RS-232 bi-directional Interface (Optional)
Calibration	Automatic External
Display	3 x 6 digits LCD digital displays
Power supply	Rechargeable battery or mains power
Operating temperature	-10°C - 40°C
Housing	ABS plastic, stainless steel pan
Overall dimensions (wxdxh)	315x355x110mm
Net weight	4.1kg

## Key Accessories

- Hard Carry Case with Lock
- In-use Cover
- In-Use Cover pack of 5
- In-Use Cover pack of 10
- Thermal Paper Roll for Adam Printer
- RS-232 Cable
- AdamDU - Data Collection Program
- Adam Printer
- RS-232 to USB Interface Cable (must be ordered with 3074010266 RS232 cable)
- ATP Thermal Printer
- Thermal Paper for ATP



Adam Equipment - offices worldwide

[www.adamequipment.com](http://www.adamequipment.com)