



Speed, Performance, Value - The Right Balance

- 3 Weighing units (g, lb, oz)
- Hold function
- Full range tare
- Simple user-friendly operation
- Large easy to read LCD display
- Simple 2 button operation
- Solid 1 piece construction is lightweight and easy to clean
- Easy grip carrying handles
- Splashproof to protect from accidental spills
- Lightweight and portable
- Low cost and simple to use
- External calibration
- Battery power (2x AA Batteries)
- Auto sleep / power down function to save battery life
- AC Adapter
- Weighing
- Animal / dynamic weighing




Technical Specifications

Maximum capacity	20kg
Readability	10g
Repeatability (S.D.)	10
Linearity (\pm)	20g

General Specifications

Model	MXB Baby
Units of measurement	Kilogram (kg), Pound (lb), Pound-Ounce (lb-oz)
Power supply	6V AC adaptor or 4 x AA Batteries
Display	LCD
Overall dimensions (wxdxh)	600x340x70mm
Net weight	2kg



Adam Equipment - offices worldwide

www.adamequipment.com