





## Speed, Performance, Value - The Right Balance

- 9 weighing units
- Date and time
- Full range tare
- Zero Tracking
- Multilingual display
- Large backlit display with dual text prompts
- Capacity tracker
- Large stainless steel pan
- Solid metal housing
- Overload protection
- Lock down mounting slot for Kensington™ type lock and cable
- Non-slip adjustable levelling feet
- Colour coded, sealed keypad
- Dual tare keys
- Removable draft shield supplied on 0.001g readabilities
- Modern low profile design
- RS-232 bi-directional interface
- GLP print outs
- External calibration
- Rechargeable battery
- AC Adapter
- Auto sleep / power down function to save battery life
- Weighing
- Parts counting
- Percentage weighing
- Check weighing
- Density determination
- Below balance weighing



## Technical Specifications

Model	PGL 203	PGL 303	PGL 2002	PGL 3002	PGL 4001	PGL 6001	PGL 8001	PGL 10001	PGL 20001
<b>Maximum capacity</b>	200g	300g	2000g	3000g	4000g	6000g	8000g	10kg	20kg
<b>Readability</b>	0.001g	0.001g	0.01g	0.01g	0.1g	0.1g	0.1g	0.1g	0.1g
<b>Repeatability (S.D.)</b>	0.002g	0.002g	0.02g	0.02g	0.2g	0.2g	0.2g	0.2g	0.2g
<b>Linearity (±)</b>	0.004g	0.004g	0.04g	0.04g	0.4g	0.4g	0.4g	0.4g	0.4g
<b>Units of measurement</b>	Gram (g), Milligram (mg), Carat (ct), Grain (GN), Newton (N), Troy Ounce (Ozt), Pennyweight (dWt), Ounce (Oz), Custom unit		Gram (g), Kilogram (kg), Carat (ct), Grain (GN), Newton (N), Troy Ounce (Ozt), Pennyweight (dWt), Ounce (Oz), Custom unit		Gram (g), Kilogram (kg), Carat (ct), Grain (GN), Newton (N), Pound (lb), Ounce (Oz), Pound-Ounce (lb-oz), Custom unit				
<b>Pan size</b>	145x125mm		192x192mm				400x300mm		
<b>Net weight</b>	4kg				6kg				

## General Specifications

<b>Typical stabilisation time</b>	2-3 Sec
<b>Operating temperature</b>	10°C 40°C
<b>Power supply</b>	AC adaptor and rechargeable battery pack
<b>Calibration</b>	External
<b>Display</b>	Backlit LCD display with dual digits and capacity tracker
<b>Overall dimensions (wxdxh)</b>	358x251x104mm

## Key Accessories

- Density Determination Kit (mg models only)
- Security Lock and Cable
- Anti-Vibration Table
- Below Balance Hanger
- In-use Cover
- Dust Cover
- Adam Printer
- AdamDU - Data Collection Program
- RS-232 Cable
- Thermal Paper Roll for Adam Printer
- F1 1kg Calibration Weight for PGL 3002
- F1 1kg Calibration Weight for PGL 2002
- F1 200g Calibration Weight for PGL 303
- F1 2kg Calibration Weight for PGL 4001
- 2 x F1 5kg Calibration Weights for PGL 20001
- F1 200g Calibration Weight for PGL 203
- F1 5kg Calibration Weight for PGL 10001
- F1 5kg Calibration Weight for PGL 8001
- F1 5kg Calibration Weight for PGL 6001



Adam Equipment - offices worldwide

[www.adamequipment.com](http://www.adamequipment.com)