

CBW Bench Check Weighing Scales



Speed, Performance, Value - The Right Balance

- Non-slip adjustable feet
- Large backlit LCD display
- Overload protection
- Easy to use keypad
- Rechargeable battery up to 70 hours operation
- Tare function
- Large stainless steel pan
- Memory accumulation
- Pre-set weighing with alarm
- Weighing
- Check weighing
- Percentage weighing
- Weight accumulation
- Parts counting



www.adamequipment.com

Technical Specifications

Model	CBW 3	CBW 6H	CBW 6	CBW 15	CBW 30	CBW 45
Maximum capacity	3000g	6000g	6000g	15kg	30kg	45kg
Readability	0.1g	0.1g	0.2g	0.5g	1g	2g
Resolution	1:30000	1:60000	1:30000	1:30000	1:30000	1:22500
Tare range	-3kg	-6kg	-6kg	-10kg	-30kg	-45kg
Repeatability (S.D.)	0.1g	0.1g	0.2g	0.5g	1g	2g
Linearity (\pm)	0.2g	0.2g	0.4g	1g	2g	4g
Pan size	225x275mm					

General Specifications

Units of measurement	Gram (g), Kilogram (kg)
Applications	Weighing Scales
Functions	Weighing, basic parts counting, % weighing and checkweighing
Typical stabilisation time	2 Sec
Calibration	Automatic External
Display	1 x 6 digits LCD digital display
Power supply	9 VDC, 800 mA from external adapter
Operating temperature	0°C - 40°C
Housing	ABS plastic, stainless steel platform
Overall dimensions (wxdxh)	315x355x110mm
Net weight	4.1kg

Key Accessories

- In-use Cover
- Adam Printer
- In-Use Cover pack of 10
- In-Use Cover pack of 5
- Weigh Below Hanger (factory fitted)
- Thermal Paper Roll for Adam Printer



Adam Equipment - offices worldwide

www.adamequipment.com