

Technical Specifications

Maximum capacity	6000g
Readability	2g
Repeatability (S.D.)	2g
Linearity (\pm)	2g
Pan size	225x275mm

General Specifications

Model	QBW 6000
Resolution	1:3000
Units of measurement	Gram (g), Kilogram (kg), Pound (lb), Ounce (Oz), Pound-Ounce (lb-oz)
Typical stabilisation time	2 Sec
Calibration	Manual external
Display	Large LCD display with 25mm/1" digits
Power supply	Rechargeable battery / mains adapter (adapter supplied as standard)
Operating temperature	0°C - 40°C
Housing	ABS Plastic, stainless steel platform
Overall dimensions (wxdxh)	315x355x110mm
Net weight	4.1kg

Key Accessories

- Weigh Below Hanger (factory fitted)
- In-use Cover
- In-Use Cover pack of 10
- Hard Carry Case with Lock
- In-Use Cover pack of 5



Adam Equipment - offices worldwide

www.adamequipment.com