

Technical Specifications


| | |
|----------------------|-----------|
| Maximum capacity | 6000g |
| Readability | 2g |
| Repeatability (S.D.) | 2g |
| Linearity (±) | 2g |
| Pan size | 225x275mm |

General Specifications

| | |
|----------------------------|--|
| Model | QBW 6000 |
| Resolution | 1:3000 |
| Units of measurement | Gram (g), Kilogram (kg), Pound (lb), Ounce (Oz), Pound-Ounce (lb-oz) |
| Applications | General purpose weighing scale |
| Typical stabilisation time | 2 Sec |
| Calibration | Manual external |
| Display | Large LCD display with 25mm/1" digits |
| Power supply | Rechargeable battery / mains adapter (adapter supplied as standard) |
| Operating temperature | 0°C - 40°C |
| Housing | ABS Plastic, stainless steel platform |
| Overall dimensions (wxdxh) | 315x355x110mm |
| Net weight | 4.1kg |

Key Accessories

- Weigh Below Hanger (factory fitted)
- In-use Cover
- In-Use Cover pack of 10
- Hard Carry Case with Lock
- In-Use Cover pack of 5



Adam Equipment - offices worldwide

www.adamequipment.com