

## Technical Specifications


Model	TBB 2610T	TBB 2610S
Maximum capacity	2610g	2610g
Readability	0.1g	0.1g
Tare range	+225g	-
Repeatability (S.D.)	0.2g	0.2g
Linearity ( $\pm$ )	0.2g	0.2g
Pan size	152mm $\varnothing$	

## General Specifications

Units of measurement	Gram (g)
Applications	Weighing, Net weight / tare, Below balance weighing, Density of liquids or solids, Filling
Calibration	External
Display	Graduated Beam
Housing	Aluminium
Overall dimensions (wxdxh)	505x110x160mm
Net weight	2.7kg

## Key Accessories

- Security Lock and Cable



Adam Equipment - offices worldwide

[www.adamequipment.com](http://www.adamequipment.com)