


Technical Specifications

Maximum capacity	311g
Readability	0.01g
Repeatability (S.D.)	0.02g
Linearity (\pm)	0.02g
Pan size	100mm \varnothing

General Specifications

Model	AQ 311S
Units of measurement	Gram (g)
Applications	Broad range
Functions	Weighing
Calibration	External
Display	Graduated Beam
Other features	Magnetic damping, Density platform with adjustable height
Housing	Aluminium
Overall dimensions (wxdxh)	370x300x137mm
Net weight	1.3kg



Adam Equipment - offices worldwide

www.adamequipment.com