


## Technical Specifications

Maximum capacity	20kg
Readability	10g
Repeatability (S.D.)	10
Linearity ( $\pm$ )	20g

## General Specifications

Model	MXB Baby
Units of measurement	Kilogram (kg), Pound (lb), Pound-Ounce (lb-oz)
Power supply	6V AC adaptor or 4 x AA Batteries
Display	LCD
Overall dimensions (wxdxh)	600x340x70mm
Net weight	2kg



Adam Equipment - offices worldwide

[www.adamequipment.com](http://www.adamequipment.com)