

## Technical Specifications

| Model                | LBK 3     | LBK 6 | LBK 12 | LBK 30 |
|----------------------|-----------|-------|--------|--------|
| Maximum capacity     | 3000g     | 6000g | 12kg   | 30kg   |
| Readability          | 0.5g      | 1g    | 2g     | 5g     |
| Repeatability (S.D.) | 0.5g      | 1g    | 2g     | 5g     |
| Linearity ( $\pm$ )  | 1g        | 2g    | 4g     | 10g    |
| Pan size             | 250x180mm |       |        |        |

## General Specifications

|                            |  |
|----------------------------|--|
| Typical stabilisation time | 2 Sec  |
| Units of measurement       | Gram (g), Kilogram (kg), Pound (lb), Ounce (Oz), Pound-Ounce (lb-oz) |
| Operating temperature      | 0° - 40°C  |
| Power supply               | 12VDC 800mA adapter and 6V 4 Ah rechargeable battery                 |
| Overall dimensions (wxdxh) | 250x307x110mm  |
| Net weight                 | 1.8kg  |

## Key Accessories

- In-use Cover
- In-Use Cover pack of 10
- In-Use Cover pack of 5



Adam Equipment - offices worldwide

[www.adamequipment.com](http://www.adamequipment.com)