



LBK *Weighing Scales*

LBK bench scales are easy to set up and use right out of the box. The sleek, low-profile design and compact size make these scales suitable for manufacturing, quality assurance and tight laboratory work spaces. Five color-coded keys highlight frequently used functions for basic weighing, percentage weighing and parts counting. LBK is an excellent value and provides a durable and practical solution for many industries.




LBK Weighing Scales

Features

- Vivid, backlit LCD easily visible in any lighting conditions
- Color-coded keys facilitate quick recognition of the most frequently used buttons
- Level indicator and adjusting feet ensure proper balance setup for optimum weighing results
- Large, type 304 stainless steel pan allows swift cleaning
- Lightweight design simplifies transport and portability
- Solid ABS housing is durable for industrial uses
- Audible overload alarm warns when capacity has been exceeded
- External calibration allows for verification and adjustment with weights
- Parts counting with preset sample sizes
- Zero-tracking feature ensures display returns to zero reading
- AC adapter included
- Rechargeable battery included for operation virtually anywhere
- Programmable auto power-off to save energy

Applications

- Weighing
- Parts counting
- Percentage weighing



Model	LBK 6a	LBK 12a	LBK 25a	LBK 65a
Capacity	6lb / 3kg	12lb / 6kg	25lb / 12kg	65lb / 30kg
Readability	0.001lb / 0.5g	0.002lb / 1g	0.005lb / 2g	0.01lb / 5g
Repeatability (S.D.)	0.001lb / 0.5g	0.002lb / 1g	0.005lb / 2g	0.01lb / 5g
Linearity (+/-)	0.002lb / 1g	0.004lb / 2g	0.01lb / 4g	0.02lb / 10g
Pan Size	9.8"x7.1" / 250x180mm			
Weighing Units	g, kg, lb, lb:oz, oz			
Stabilization Time (sec)	2			
Calibration	External calibration			
Display	Backlit LCD with 0.9" / 24mm-high digits			
Power Supply	12VDC 800mA adapter / Internal rechargeable battery			
Operating Temperature	0° to 40°C / 32° to 104°F			
Housing	ABS			
Overall Dim.	9.8"x12.1"x4.3" / 250x307x110mm			
Net Weight	3.97lb / 1.80kg			

For more than 40 years, professionals worldwide have relied on Adam Equipment for an extensive selection of dependable and affordable weighing equipment. Adam provides the right balance of speed, performance and value.

ACCESSORIES

- | | |
|-----------|-------------------------------|
| 700200056 | In-use wet cover |
| 700200063 | In-use wet cover (pack of 5) |
| 700200064 | In-use wet cover (pack of 10) |

Supplied by: